

IV-161 DIGITAL CONDUCTIVITY METERS

PRODUCT FEATURES

- Inbuilt stirrer
- Bright LED display
- Stable and accurate result
- Easy to operate

ACCESSORIES

- Magnetic capsule
- Conductivity cell 1.0
- Rotational conductivity cell stand
- Instruction manual
- Dust cover

INBUILT STIRRER



Technical specification

Display	3½ Digit LED
Measurement	Conductivity
Cond. Ranges	Five Ranges
Range Selection	Automatic 5 Ranges
Cond.	±0.5% FS ±1DIGIT
Temperature Compensation	Manual 0 to 50°C
Stirrer	Inbuilt Stirrer (With Speed Variable)
Cell Constant	Adjustable on Digital Display
Measuring Cell	Platinum DIP Type
Cell Const. Measurement	Auto
Power	230V ± 10% AC
Accessories	Conductivity Cell, Magnetic Capsule, Operation Manual. Dust Cover, Cell Stand