

## IV-12 DIGITAL PHOTO COLORIMETER

### PRODUCT FEATURES

- Bright LED display
- Accurate and stable result
- Set zero with potentiometer adjustment
- Long life silicon photodiode

### ACCESSORIES

- Glass test tubes - 5 Nos.
- Instruments cover - 1 Nos.
- Operating Manual - 1 Nos.
- Test tube stand - 1 Nos.
- Dust cover - 1 Nos.



### Technical specification

<b>Filters</b>	8 Filters
<b>Accuracy</b>	±0.02 O.D.
<b>Reproducibility</b>	±0.01 O.D.
<b>Stability</b>	±0.02 O.D. Per Hour
<b>Light Source</b>	3.5V LED Bulb
<b>Key Board</b>	Potentiometer Adjustment
<b>Amplifier</b>	I.C. Based
<b>Detector</b>	Selenium Photocell/photodiode
<b>Warm-up Time</b>	5 Minutes
<b>Sample System Tubes</b>	12mm Dia Matched Glass Test
<b>Power Source</b>	230V±10% 50hz.ac
<b>Dimension</b>	280x260x160mm approx
<b>Weight</b>	3kg