


## IV-14 PORTABLE INSTRUMENTS PH METER



### Technical specification

<b>Model No.</b>	
<b>Model No.</b>	IV-14
<b>Range</b>	0-14 pH
<b>Resolution</b>	0.01 pH
<b>Impedance</b>	$10^{12}$ ohms
<b>Sensor</b>	combination pH electrode
<b>Battery</b>	9V DC
<b>Accuracy</b>	0.01 pH
<b>Temp. Comp.</b>	0-100° C (Manual)
<b>Display</b>	3½ digit LCD
<b>Size</b>	180 x 85x 30 mm
<b>Weight</b>	250 gm (Approx)