

## LT-92 OVER HEAD STIRRERS


Overhead Stirrer is suitable for mixing and stirring low, medium and high viscosity liquids in the laboratory. The stirring machine adopts high-quality engineering plastic as the body, the operating is stable, the output uses variable speed increase mechanism, so that the output torque is multiplied.

### PRODUCT FEATURES

- The LCD screen displays the set speed and actual speed, which is intuitive and clear
- DC permanent magnet motor, low noise, maintenance-free, stable operation
- The overload will stop automatically to ensure safe operation
- Stable start-up, effectively preventing sample overflow
- Automatically memorize the final set working parameters, multi-segment programming control



### Technical specification

Image	
Product Model	LT-92
Voltage V	AC100-240
Frequency Hz	50-60
Power w	72
Timing range	1-99h 59 min
Display method	LCD
Speed Range rpm	100-2000
Speed Control Accuracy rpm	±5%
Max Stirring Volume (H2O) L	20
Recommended Mixing Amount (H2O) L	5
Stirring Max Viscosity mPa s	10000
Max Torque N. cm	10
Motor type	DC Permanent Magnet Motor