

IV-47 MINI CENTRIFUGE


Mini Centrifuge equipped with a variety of centrifuge rotors and a variety of test tube sets. With flip open function, automatically stop when the cover is opened. It can be used for serum extraction from whole blood, various sample extraction, rapid sample sedimentation, separation of trace blood cells, microbial sample processing, PCR experiment zone centrifugation, etc.

PRODUCT FEATURES

- The unique rotor can meet different experiments without change.
- Micro DC motor and carbon fiber rotor.
- With silent technology operate with ultra-low noise.
- Strong shell made of high-strength transparent material, resistance to fall and stock.
- Flip switch function, automatically stop when the cover is opened.



Technical specification

Image	
Product Model	IV-47
Digital Display	Yes
Speed Range RPM	6000
Controller Type	Adjustable speed & Timing
Size	180x118x160mm
Power Consumption	9w
Input Voltage	24vdc 1.5 amp
Motor type	Dc Permanent Magnet Motor
Weight	500 gm